

# PGSE Mini Servo Gripper User Guide

---

Issue: V1.0  
Date: 2024-08-22

**Copyright © Shenzhen Yuejiang Technology Co., Ltd 2024. All rights reserved.**

No part of this document may be reproduced or transmitted in any form or by any means without the prior written consent of Shenzhen Yuejiang Technology Co., Ltd (hereinafter referred to as “Dobot”).

### **Disclaimer**

To the maximum extent permitted by applicable law, the products described (including its hardware, software and firmware, etc.) in this document are provided **AS IS**, which may have flaws, errors or faults. Dobot makes no warranties of any kind, express or implied, including but not limited to, merchantability, satisfaction of quality, fitness for a particular purpose and non-infringement of third party rights. In no event will Dobot be liable for any special, incidental, consequential or indirect damages resulting from the use of our products and documents.

Before using our product, please thoroughly read and understand the contents of this document and related technical documents that are published online, to ensure that the robot arm is used on the premise of fully understanding the robot arm and related knowledge. Please use this document with technical guidance from professionals. Even if follow this document or any other related instructions, damages or losses may happen in the using process. Dobot shall not be considered as a guarantee regarding all security information contained in this document.

The user has the responsibility to make sure following the relevant practical laws and regulations of the country, in order that there is no significant danger in the use of the robot arm.

## **Shenzhen Yuejiang Technology Co., Ltd.**

Address: Room 1003, Building 2, Chongwen Garden, Nanshan iPark, Liuxian Blvd, Nanshan District, Shenzhen, Guangdong Province, China

Website: [www.dobot-robots.com](http://www.dobot-robots.com)

## Preface

### Purpose

This document introduces the functions, parameters, installation, and operation of the PGSE Mini Servo Gripper, helping users understand and use the servo gripper effectively.

### Intended audience

This document is intended for:





- Customer
- Sales Engineer
- Installation and Commissioning Engineer
- Technical Support Engineer

### Revision history

Date	Version	Revised content
2024/08/22	V1.0	The first release

### Symbol conventions

The symbols that may be found in this document are defined as follows.

Symbol	Description
 DANGER	Indicates a hazard with a high level of risk which, if not avoided, could result in death or serious injury.
 WARNING	Indicates a hazard with a medium level or low level of risk which, if not avoided, could result in minor or moderate injury, robot arm damage.
 NOTICE	Indicates a potentially hazardous situation which, if not avoided, could result in robot arm damage, data loss, or unanticipated result.
 NOTE	Provides additional information to emphasize or supplement important points in the main text.

## Contents

<b>Preface.....</b>	<b>ii</b>
<b>1. Product Introduction .....</b>	<b>1</b>
1.1 Product specifications .....	1
1.2 Indicator definition.....	1
1.3 Mechanical specifications.....	2
1.4 Electrical interface .....	2
<b>2. Installation Instructions .....</b>	<b>3</b>
2.1 Assembling the gripper .....	3
2.2 Installing gripper onto Magician E6 end.....	4
<b>3. Gripper Control .....</b>	<b>5</b>

## 1. Product Introduction

### 1.1 Product specifications

The PGSE157 is a parallel servo gripper equipped with a pair of parallel fingertips that move symmetrically during operation. The gripper comes with an aviation connector cable that can be directly connected to the Magician E6 end. It communicates with Magician E6 via digital I/O.



Parameter	Value
Model	DT-AC-PGSE157-01
Gripping force (per jaw)	6 N – 15 N
Stroke	7 mm
Maximum recommended payload*	0.25 kg
Opening/Closing time	0.15 s / 0.15 s
Size	85.6 mm x 38 mm x 23.2 mm
Weight	0.15 kg
Noise emission	< 50 dB
Protection level	IP40
Communication protocol	Digital I/O
Operating voltage	24V DC $\pm$ 10%
Rated current	0.15A
Peak current	0.8A
Rated power	3.6W
Recommended environment	0°C – 40°C, $\leq$ 85% RH

\* The maximum recommended payload depends on the shape of the object being gripped, the material of the contact surface, friction, and motion acceleration. The center of gravity of the gripped object can also affect the payload. Please consult us if you have any questions.


### 1.2 Indicator definition

The gripper provides real-time feedback on its status. In addition to reading through digital I/O, the status can also be judged by the indicator light color, described as follows:

- Uninitialized status: Red light flashing, other lights OFF.
- Initialized status: Blue light ON, indicating the gripper is in operational status.


- Object gripped status: Green light ON, other lights OFF.

### 1.3 Mechanical specifications



### 1.4 Electrical interface

The gripper connects to the end of Magician E6 through an aviation connector, with DI/DO in PNP type. The interfaces are defined as follows:




Pin	Definition
1	24V
2	DO1
3	DO2
4	DI1
5	DI2
6	GND


## 2. Installation Instructions

### 2.1 Assembling the gripper

1. Attach the L-shaped gripper adapter to the electromagnetic gripper and tighten the screws.




2. Install the gripper fingers and tighten the screws. Adjust the installation position of the gripper fingers on the L-shaped adapter according to the actual situation. The protective pad on the side of the gripper fingers can be attached or removed as needed.



## 2.2 Installing gripper onto Magician E6 end

1. Attach the Magician E6 adapter flange (an accessory of Magician E6) to the gripper and tighten the two screws.



2. Install the electromagnetic gripper onto the Magician E6 end through the adapter flange and tighten the four screws.



3. Connect the gripper cable to the end-of-arm aviation connector.



### 3. Gripper Control

The robot can control the gripper through Tool DO1 and DO2.

The gripper has been initialized at the factory. In most cases, users can directly control the gripper using Tool DO1.

Tool_DO1	Function	Tool_DO2	Function
OFF	Open	OFF	Default
ON	Close	ON	Initialize

The robot can read the gripper's status through Tool DI1 and DI2.

Tool_DI1	Tool_DI2	Gripper status
OFF	OFF	Moving/Just started/Fault
ON	OFF	Closed
OFF	ON	Opened
ON	ON	Object gripped