


Topgene & Osmunda Medical Device Co., Ltd.

Introduction & Documents

Surgical
FACE MASK



Environment

Factory Appearance




Production Line



Product Specification

- Mask



- Packaging



- **Product Details**

Product name: Surgical Face Mask

Model: Elastic Earloop (175*95mm)

Materials: The face mask is composed of 3 ply fabric. Outer layer and inner layer are made of SSS non-woven fabric and the middle layer with a filter paper.

Application range: Protection for medical organization in the general medical environment

Performance: BFE, Delta P conform to the standards of EN14683:2019 and YY/T 0969-2013

Package: 50pcs per box / 2000pcs per container

Certification: NMPA (China), CE (Europe)

European reference standard: Type IIR

- **Instructions**

1. Take out the mask and apply it on your face directly. Pull open the folded part until it completely covers the bridge of nose and the chin.
2. Press down both sides of the nose bridge bar to match the nose bridge, and adjust the mask belt to fit the face shape.

- **Notes**

1. Please check the packing before use, And use it within the validity period.
2. This product is for single use. Do not reuse it.
3. The person who is non-woven fabric allergy should use with caution.

- **Storage**

Store in a dry, ventilated environment with no corrosive gases.

- **Validity period: 2 years**